Part1:- Read in and subsample CyTOF data

Open the fcs files within Python and merges the files in a single step w/o having to open the files in Fcs Express or Flowjo.

Import packages and enter path to directory with all the data

```
In [ ]: | from FlowCytometryTools import FCMeasurement
        import os
        import pandas as pd
        pd.set_option("max_columns", 50)
        path = input(r"Enter path: ")
        os.chdir(path)
```

Open files and check if multiple panels are available for analysis

```
In [3]: panels = []
           for f in os.listdir():
    if ".fcs" in f:
                     panels.append(f.split('_')[3])
           panels = list(set(panels))
           print("Panels available:{}" . format(panels))
term = input(r"Select panel to proceed: ")
           Panels available:['Tcell']
           Select panel to proceed: Tcell
```


Open files of the selected panel and subsample the events

Subampling done such that all samples have the same number of events, n = smallest sample size by defaullt

```
In [5]: sample_size = []
          num_files = 0
          for f in os.listdir():
               if term in f and f.endswith('.fcs'):
                    #print(f)
                   datafile = f
sample = FCMeasurement(ID = 'sample1', datafile= datafile)
                   sample_size.append(len(sample.data))
                   num_files += 1
          print("{} files found".format(num_files))
print("minimum sample size is {}".format(min(sample_size)))
          5 files found
```

minimum sample size is 2163

```
In []: # This will need editing based on the names of the columns and file naming scheme to get tissue name
          answer = input(r"Type Y to proceed or N to exit: ")
          if answer == 'Y':
               print("proceeding to subset the data to include {} cells from each site".format(min(sample_size)))
               opts = input(r"Change default sample size (Y/N): ")
                    merged_df = pd.DataFrame()
                    for f in os.listdir():
                         if term in f and f.endswith('.fcs'):
                              datafile = f
                              sample = FCMeasurement(ID = 'sample1', datafile= datafile)
                              sample_size.append(len(sample.data))
                              parameter = sample.channels['$PnS']
                              mod_parameters = []
data_columns = []
                              # Here we are cleaning up the column names (markers) to generate user friendly names
# Edit this based on how column names are defined in your sample
                              # Lot this based on now cotamn names are defined in your sample
for p in parameter:
    if '_' in p and 'dist' not in p and 'Time' not in p and 'Event' not in p and 'DNA' not in p \
    and 'Osmium' not in p and 'Viability' not in p and 'PD_1' not in p:
        mod_parameters.append(p.split('_', maxsplit = 1)[1])
        if len(p.split('_', maxsplit = 1)) > 1 and 'PD_1' not in p:
                                   data_columns.append(p.split('_', maxsplit = 1)[0])
elif 'PD_1' in p:
                                        data_columns.append(p)
                                        mod_parameters.append(p)
                                   elif '_ not in p and 'dist' not in p and 'Time' not in p and 'Event' not in p and 'DNA' not in
          p \
                                   and 'Osmium' not in p and 'Viability' not in p:
                                        data_columns.append(p)
                                        mod_parameters.append(p)
                                   else:
                                        mod_parameters.append(p)
                              sample_data = pd.DataFrame(sample.data)
                              sample data.columns = mod parameters
                              sub_data = sample_data.sample(n = min(sample_size))
# Get tissue name, edit based on file naming scheme
sub_data['tissue'] = f.split('_')[4]
                              merged_df = pd.concat([merged_df, sub_data], axis = 0)
                    'GranzymeB']
                    merged_df.reset_index(inplace = True, drop = True)
                    print("data merged from {} files".format(num_files))
                    print("merged dataset size has {} events".format(len(merged_df)))
                     #merged_df.head()
               if opts == 'Y':
                    .
sub_num = input(r"Enter number of events to subsample: ")
                    sub_num = int(sub_num)
                    merged_df = pd.DataFrame()
                    for f in os.listdir():
                         if term in f:
                              sample = FCMeasurement(ID = 'sample1', datafile= datafile)
                              sample_size.append(len(sample.data))
                              parameter = sample.channels['$PnS']
                              mod_parameters = []
                              data_columns = []
                              # Here we are cleaning up the column names (markers) to generate user friendly names
                              # Edit this based on how column names are defined in your sample
                              for p in parameter:
                                   p in parameter.
if '_' in p and 'dist' not in p and 'Time' not in p and 'Event' not in p and 'DNA' not in p \
and 'Osmium' not in p and 'Viability' not in p and 'PD_1' not in p:
    mod_parameters.append(p.split('_', maxsplit = 1)[1])
    if len(p.split('_', maxsplit = 1)) > 1 and 'PD_1'not in p:
                                             data_columns.append(p.split('_', maxsplit = 1)[0])
                                   elif 'PD_1' in p:
                                        data_columns.append(p)
                                   mod_parameters.append(p)
elif '_' not in p and 'dist' not in p and 'Time' not in p and 'Event' not in p and 'DNA' not in
          p \
                                   and 'Osmium' not in p and 'Viability' not in p:
                                        data_columns.append(p)
                                        mod_parameters.append(p)
                                   else:
                                        mod_parameters.append(p)
                              sample_data = pd.DataFrame(sample.data)
                              sample_data.columns = mod_parameters
                              sub_data = sample_data.sample(n = sub_num)
                              # Get tissue name. edit based on file naming scheme
```



Part 2:- Data transformation and dimensionality reduction

Arcsin transform the data and run dimensionality reduction using TSNE and UMAP. Then run Louvain clustering to identify cluster of cells

```
In [66]: import matplotlib.pyplot as plt
         import numpy as np
         import phenograph
         import umap
         import seaborn as sns
         import matplotlib.patches as mpatches
         from matplotlib import colors as mcolors
         from MulticoreTSNE import MulticoreTSNE as TSNE
         from scipy import stats
In [67]: # Transforming the data
         cof = 5 # cofactor for arcsin transformation
         trans_data = np.arcsinh(merged_df[sel_cols]/cof)
         trans data.head()
Out[67]:
               Thet
                     TIGIT
                            HI ADR
                                    CD123
                                            CD27
                                                    Bcl6
                                                           CD1c CD45RA
                                                                          PD 1
                                                                                 CD69
                                                                                        CD56
                                                                                               CCR7
                                                                                                      CD103
                                                                                                             CXCR5
                                                                                                                      41RR
                                                                                                                            CX
         0 0.000000 0.000000 1.696996 0.000000 0.565838 4.548936 0.000000
                                                                1.132587 0.109828 0.224523 0.153296
                                                                                                                           0.25
                                                                                             0.000000
                                                                                                    0.437203
                                                                                                            0.099043
                                                                                                                    0.000000
         1 1.511604 0.288182 1.172662 0.173042 0.100105 1.709869 0.671835 2.893921 0.000000 2.505974 0.298781 0.140513 0.598210 0.017804
         2 1.161298 0.551242 0.000000 0.650896 0.000000 3.689836 0.461862 0.000000 0.000000 2.944076 2.505732 0.000000 3.613437 0.225402 0.000000 1.87
          4 0.000000 0.000000 1.002638 1.031823 0.193457 4.462292 0.000000 0.000000 0.122732 0.435399 0.000000 0.184861 0.329938 0.000000 0.000000 0.77
```

Dimensionality reduction using TSNE

```
In [68]: | tsne = TSNE(n_components=2, perplexity=30.0, early_exaggeration=12.0, learning_rate=200.0,
                      n_iter=1000, n_iter_without_progress=300, min_grad_norm=1e-07, metric='euclidean'
                      init='random', verbose=0, random_state=None, method='barnes_hut', angle=0.5, n_jobs=4)
          tsne_embedding = tsne.fit_transform(trans_data.values)
In [69]:
         t_data = pd.DataFrame(tsne_embedding, columns = ['t_x', 't_y'])
          t data.head()
Out[69]:
                 t x
                          t_y
             3.597770 -25.014189
             9.390791
                      5.471216
            22.845692 24.592626
          3 15.841476 -15.495996
          4 -1.276508 -26.521527
```

Dimensionality reduction using UMAP

Run phenograph clustering on data

```
In [72]: | communities, graph, Q = phenograph.cluster(trans_data.values, k = 150)
           num_clusters = len(set(communities)) - (1 if -1 \frac{1}{1} communities else 0)
          print(num_clusters)
          Finding 150 nearest neighbors using minkowski metric and 'auto' algorithm
          Neighbors computed in 7.4553444385528564 seconds
          Jaccard graph constructed in 122.9233911037445 seconds
          Wrote graph to binary file in 17.133769750595093 seconds Running Louvain modularity optimization
          After 1 runs, maximum modularity is Q = 0.817754
          After 3 runs, maximum modularity is Q = 0.818823
Louvain completed 23 runs in 27.582430124282837 seconds
          PhenoGraph complete in 175.63334846496582 seconds
In [73]: post_df = pd.concat([t_data, umap_df], axis=1)
           post_df['cluster'] = communities
           post_df.head()
Out[73]:
                                                u_y cluster
                              t_y
           0 3.597770 -25.014189
                                   5.866962 4.581736
                                                         2
              9.390791
                         5.471216
                                  -7.496478 1.218129
                                                        10
           2 22.845692 24.592626 -10.545886 3.798560
                                                         9
                                                         6
           3 15.841476 -15.495996
                                   5.594634 2.337136
           4 -1.276508 -26.521527
                                   5.804407 5.434681
                                                         0
```

Visualize the clusters

Define colors that will be used to color clusters

```
In [74]: colors = dict(mcolors.BASE_COLORS, **mcolors.CSS4_COLORS)
palette = []
for k, v in colors.items():
    if '#' in v:
    palette.append(v)

palette = palette[0::4]
del(palette[0])
del(palette[0])

In [75]: clust_color = dict(zip(sorted(post_df['cluster'].unique()), palette))
sample_color = post_df['cluster'].map(clust_color)
```


```
In [76]: # Add color and tissue information to post dim-reduction dataframe
          post_df['color'] = sample_color
          post_df['Tissue'] = merged_df['tissue']
          post_df.head()
Out[76]:
                                    u_x
                                            u_y cluster
                                                          color Tissue
                           t_y
          0 3.597770 -25.014189
                                5.866962 4.581736
                                                    2 #FF7F50
                                                               LUNG
          1 9.390791 5.471216 -7.496478 1.218129
                                                   10 #FF00FF LUNG
          2 22.845692 24.592626 -10.545886 3.798560
                                                    9 #1E90FF LUNG
          3 15.841476 -15.495996 5.594634 2.337136
                                                    6 #8B0000 LUNG
          4 -1.276508 -26.521527 5.804407 5.434681
                                                    0 #FFEBCD LUNG
```

Generate the scatter plot for all samples


```
In [77]: def generate_all_scatter(df, embedding):
                     if embedding == 'tsne':
                           dr = 't
                     if embedding == 'umap':
                           dr = 'u
                     # Define x and y variables for the plot
                     x = df['{}_x'.format(dr)].values
y = df['{}_y'.format(dr)].values
                     clr = df['color']
                     sns.set_style('whitegrid')
                     sns.set_styte( whitegrid )
plot_kwds = {'alpha' : 0.4, 's' : 15, 'linewidths':0}
figure, ax = plt.subplots(figsize=(10,9))
ax.set_facecolor('k')
ax.scatter(x, y, c= clr,**plot_kwds)
clust_col_new = {i:clust_color [i] for i in sorted(df.cluster.unique())}
                     # Add labels to each cluster
                     for i, label in enumerate(clust_col_new.keys()):
                           plt.annotate(label,
                                                df.loc[df['cluster']==label,['{}_x'.format(dr),'{}_y'.format(dr)]].mean(),
horizontalalignment='center', verticalalignment='center',
size=18, weight='bold', color = 'white')
                     frame = plt.gca()
                     frame = ptt.gea()
frame.axes.get_xaxis().set_visible(False)
frame.axes.get_yaxis().set_visible(False)
                     plt.title('{}\nEstimated number of clusters: {}'.format(embedding, num_clusters),fontsize=20)
```

```
In [78]: # Call function, provide df with dimensionality reduction results and embedding to plot
generate_all_scatter(post_df, 'tsne')
```

tsne
Estimated number of clusters: 13




Generate tissue-wise scatter plots




Frequency of each cluster in tissues

```
In [83]: # Expression data for all the cells
        ex_data = pd.concat([trans_data, post_df], axis =1)
dict_list = []
            for site in tissues:
                t df = df[df['Tissue'] == site]
                obs = t_df.cluster.value_counts()
                 freq_dict = dict(zip(obs.index, obs.values))
            dict_list.append(freq_dict)
obs_df = pd.DataFrame(dict_list).T
            obs_df.columns = tissues
             # number of observations per site
            total_obs = obs_df.sum(axis=0, skipna=True).tolist()
             # calculate percent of each cluster per site
            for site, t_obs in zip(tissues, total_obs):
                obs_df[site] = obs_df[site].apply(lambda x: (x/t_obs)*100)
            obs_df.fillna(0, inplace = True)
            # colors for each cluster
clust_color = dict(zip(sorted(df['cluster'].unique()), palette)) # use palatte that we defined previously
            colors = obs_df.index.map(clust_color)
            # generate a cluster map for the data
            sns.set(font_scale=1.2)
            col_colors=colors)
            g.cax.set_visible(False)
            \verb"plt.setp"($\overline{g}$.ax\_heatmap.yaxis.get_majorticklabels(), rotation=0)
```



Bar plot for marker expression levels per cluster



Part3:- Identification of cell clusters of interest


Generate Heatmap - first heatmap to select cell clusters of interest

Average expression dataframe

```
In [90]: # Average expression
    avg_exp = ex_data.groupby('cluster').mean()
    avg_exp['cluster'] = avg_exp.index
    avg_exp['cluster color'] = avg_exp['cluster'].map(clust_color)

In [91]: # Write the expression dataframe with cluster info, dim_resuction, tissue info and cluster color to disk
    ex_data['Tissue'] = merged_df['tissue']
    ex_data.to_csv("Expression_data_dim_clust_tissue.csv", index = False, sep = ',')
```

Plot heatmap






Generate density plots for markers in clusters - may also help identify cells of interest

Define function to generate ridge plots

```
In [94]: import bokeh
           import colorcet as cc
           from numpy import linspace
           from scipy.stats.kde import gaussian_kde
           from bokeh.io import output_notebook, output_file, show
           from bokeh.models import ColumnDataSource, FixedTicker, PrintfTickFormatter
           from bokeh.plotting import figure
from bokeh.sampledata.perceptions import probly
           bokeh.io.reset_output()
#output_file()
           output_notebook()
           # Function to generate the plots
           def ridge(category, data, scale=3):
                return list(zip([category]*len(data), scale*data))
           def plot_ridge(df, marker):
    if marker not in ex_data.columns.tolist():
                     print("This marker is not available in dataset!\nPlease select another marker.")
                else:
                     clusters = sorted(df['cluster'].unique())
                     cluster_names = ["Cluster number " + str(c) for c in sorted(df['cluster'].unique())]
                     x = linspace((df[marker].values.min() - 20),(df[marker].values.max() + 20), 500)
                     source = ColumnDataSource(data=dict(x=x))
                     p = figure(y_range=cluster_names, plot_width=900, plot_height = 600,
    x_range=(-2, (df[marker].values.max() + 2)),
                                   toolbar_location="above", title = "{} expression across clusters".format(marker))
                     for c_num, c_name in zip(clusters, cluster_names):
                          data = df[df['cluster'] == c_num]
data = data[marker].values
                          pdf = gaussian_kde(data)
                          y = ridge(c_name, pdf(x))
                          source.add(ȳ, c_name)
p.patch('x', c_name, color = palette[c_num], alpha=0.6, line_color="black", source=source)
                     p.outline_line_color = None
                     p.background_fill_color = "#efefef"
                     p.ygrid.grid_line_color = None
p.xgrid.grid_line_color = "#dddddd"
p.xgrid.ticker = p.xaxis.ticker
                     p.axis.minor_tick_line_color = None
p.axis.major_tick_line_color = None
p.axis.axis_line_color = None
                     p.y_range.range_padding = 0.2
                     show(p)
```

(http:BökehdS.2:0)1 successfully loaded.




Part 4:- Subset data based on clusters of interest

Enter list of clusters of interest, can enter list of list of clusters if you want to suset two cell types


Scatter plots


Plot 2 varibles at a time - enter values manually

In [117]: # Call function and provide x-axis marker and y-axis marker
twoD_plot(cell_dict['Bcells'], 'HLADR', 'CD69')



In [109]: twoD_plot(cell_dict['Tcells'], 'CD4', "CD8")





Plot 2 variable at a time from a list of markers (arrange in a grid)

```
In [112]: def auto_twoD_plots(df, markers):
                                                   from itertools import combinations, permutations
                                                   marker_combs = list(combinations(markers, 2))
                                                  # Need to run this to generate grid row and col location rows = list(range(0, (len(markers)-1))) col = \theta
                                                   sub_rows = []
                                                   sub cols = []
                                                   for row in rows:
                                                                 while col < len(rows) - row:</pre>
                                                                               sub_rows.append(row)
                                                                               sub_cols.append(col)
                                                                               col +=1
                                                                 col = 0
                                                   # Generate the plots
                                                   fig, \ ax = plt.subplots((len(markers)-1), (len(markers)-1), \ figsize = (((len(markers)-1)*4), ((len(markers)-1)*4), (len(markers)-1)*4), (len(markers)-1
                                                   plt.subplots_adjust(wspace=0.3, hspace=0.3)
                                                   for row, col, comb in zip(sub_rows, sub_cols, marker_combs):
                                                                x = df[comb[0]]
                                                                y = df[comb[1]]
c = df['color']
                                                                 ax[row,col].scatter(x = x,y = y,c = c, s = 8)
ax[row,col].set_facecolor('black')
                                                                 ax[row,col].set_xlabel(comb[0], size = 12)
                                                                 ax[row,col].set_xlim(0,max(x))
                                                                 ax[row,col].set_ylabel(comb[1], size = 12)
                                                                 ax[row,col].set_ylim(0,max(y))
                                                   plt.tight_layout()
                                                   plt.show()
```

